

SureSigns Patient Monitoring Portfolio

Product comparison



VM Series	VM4 – 863063	VM6 – 863064	VM6 – 863065	VM8 – 863066	VM8 – 863068
NBP, SpO ₂	✓	✓	✓	✓	✓
ECG	3-lead	3- or 5-lead selectable	3- or 5-lead selectable	3- or 5-lead selectable	3- or 5-lead selectable
Respiration	✗	✓	✓	✓	✓
Basic arrhythmia	✗	✓	✓	✓	✓
IBP (one channel)	✗	✗	✓	✓	✓
Predictive temperature	●	✗	✗	✗	✗
Continuous temperature**	✗	✓	✓	✓	✓
Microstream CO ₂	✗	✗	✗	✓	✗
Display	8.4" SVGA TFT	8.4" SVGA TFT	8.4" SVGA TFT	10.4" SVGA TFT	10.4" SVGA TFT
Traces (max waveforms)	2	4	4	4	4
Bed rail hanger	✗	✓	✓	✓	✓
Recorder	●	●	●	✓	●

✓ Standard

● Optional

✗ Not available

** Quick (predictive) / oral temperature

PHILIPS



C1 Series	863051	863052	863053	863054
SpO ₂	✓	✓	✓	✓
NBP	✓	✓	✓	✓
Temperature	✗	✓**	✗	✓**
ECG	✗	✗	✓	✓
Recorder	● external	● external	● external	● external

** Quick (predictive) / oral temperature



MP20 Junior	M8001A #M20
SpO ₂	✓
NBP	✓
Temperature	● up to 2
ECG	✓
IBP	● up to 2
CO ₂	●
Recorder	●



VS1 Series	863055	863056	863057	863058	863059	863060	863061	863062
NBP	✓	✓	✓	✓	✓	✓	✓	✓
Temperature	✗	✓**	✗	✓**	✗	✓**	✗	✓**
SpO ₂	✗	✗	✗	✗	✓	✓	✓	✓
Recorder	✗	✗	✓	✓	✓	✓	✗	✗

** Quick (predictive) / oral temperature

✓ Standard ● Optional ✗ Not available



© 2006 Koninklijke Philips Electronics N.V.
 All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.
 Philips Medical Systems Nederland B.V. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. Microstream is a registered trademark of Oridion Medical, LTD.

Philips Medical Systems is part of Royal Philips Electronics

www.medical.philips.com
 medical@philips.com
 fax: +31 40 27 64 887

Printed in The Netherlands
 4522 962 09451/862 * JAN 2006

Philips Medical Systems
 Global Information Center
 P.O. Box 1286
 5602 BG Eindhoven
 The Netherlands